

Intact™ Specimen Capture Volumes

Choose the Best Sample Capture Size

The *Intact* Breast Lesion Excision System (*Intact* BLES) offers four different capture basket sizes. Each deploys from the same 6mm diameter wand.



<i>Intact</i> Wand diameter:	10mm	12mm	15mm	20mm
Basket dimensions (nominal):	10mm x 14mm	12mm x 17mm	15mm x 21mm	20mm x 25mm
Specimen weight:	.8g	1.1g	2.1g	3.0g

Compare *Intact* BLES Specimens* to Core Needles...

11 gauge Vacuum-Assist

8 gauge Vacuum-Assist

Specimen actual size:



Specimen dimensions:	2.4mm x 19mm	3.3mm x 22mm
Specimen weight:	0.09 grams	0.19 grams

...in Total Specimen Volume...

Total Volume Equivalent: (Number of Core Strips)	<i>Intact</i> BLES 10mm:	<i>Intact</i> BLES 12mm:	<i>Intact</i> BLES 15mm:	<i>Intact</i> BLES 20mm:
8 gauge Vacuum-Assist:	2.9×	5.0×	9.9×	20.9×
11 gauge Vacuum-Assist:	6.4×	10.9×	21.6×	45.7×

...and in Target Volume

Target Volume Equivalent: (Number of Core Strips)	<i>Intact</i> BLES 10mm:	<i>Intact</i> BLES 12mm:	<i>Intact</i> BLES 15mm:	<i>Intact</i> BLES 20mm:
8 gauge Vacuum-Assist:	9×	13×	21×	37×
11 gauge Vacuum-Assist:	17×	25×	39×	69×

*Specimen weights are based on the assumption that 1cc of breast tissue weighs about 1 gram.